

Anti-Fah Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P16930

Item Number: CM0010552

Introduction

Rabbit monoclonal antibody specific for human FAH (fumarylacetoacetase), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for studying tyrosine metabolism and liver/kidney disorders.

[Learn More](#)

Parameter	Value
Product Name	FAH Rabbit mAb
Species	Human
GeneID	2184
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-419 of human FAH (NP_000128.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16930
Protein Weight	46kDa
Shipping	Ice bag