

Anti-Fah Rabbit Polyclonal Antibody For Wb, Elisa - P16930

Item Number: CM0015210

Introduction

FAH Rabbit pAb is a polyclonal antibody raised in rabbit against full-length human FAH. Validated for WB and ELISA, it cross-reacts with human, mouse, and rat FAH. Ideal for studying tyrosine catabolism and liver diseases.

[Learn More](#)

Parameter	Value
Product Name	FAH Rabbit pAb
Remarks/Alias	FAH
Species	Human
GeneID (Human)	2184
GeneID	2184
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-419 of human FAH (NP_000128.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16930
Protein Weight	46kDa
Shipping	Ice bag