

## Anti-Faim3 Rabbit Polyclonal Antibody For Wb, Elisa - O60667

Item Number: CM0030545



### Introduction

Affinity-purified rabbit polyclonal antibody against human FAIM3 (TOSO, FcμR), validated for Western blot (WB) and ELISA. Cross-reacts with mouse and rat; target molecular weight 43 kDa. Ideal for studying IgM receptor biology in B cells.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | FAIM3 Rabbit pAb                                                                                                  |
| Remarks/Alias    | TOSO; FAIM3; FcμR                                                                                                 |
| Species          | Human                                                                                                             |
| Gene ID (Human)  | 9214                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 18-251 of human FAIM3 (NP_005440.1) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | O60667                                                                                                            |
| Protein Weight   | 43kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |