

Anti-Fam107A Polyclonal Antibody For Wb, Elisa - 095990

Item Number: CM0020310

Introduction

Rabbit polyclonal antibody against human FAM107A (DRR1/TU3A), validated for WB and ELISA; cross-reacts with mouse and rat. Ideal for studying actin dynamics, synaptic function, and cancer.

[Learn More](#)

Parameter	Value
Product Name	FAM107A Rabbit pAb
Remarks/Alias	DRR1; TU3A; FAM107A
Species	Human
GeneID (Human)	11170
GeneID	11170
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-144 of human FAM107A (NP_009108.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O95990
Protein Weight	17kDa
Shipping	Ice bag