

Anti-Fam107A Rabbit Polyclonal Antibody For Wb, Elisa - O95990

Item Number: CM0019921

Introduction

Rabbit polyclonal antibody targeting human FAM107A. Suitable for Western blot and ELISA; cross-reacts with human and mouse. Ideal for studying actin dynamics, synaptic function, and tumor-related pathways.

[Learn More](#)

Parameter	Value
Product Name	FAM107A Rabbit pAb
Remarks/Alias	DRR1; TU3A; FAM107A
Species	Human
Gene ID	11170
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-144 of human FAM107A (NP_009108.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O95990
Protein Weight	17kDa
Shipping	Ice bag