

Anti-Fam117A Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9C073

Item Number: CM0020526

Introduction

This rabbit polyclonal antibody targets FAM117A, a C/EBP-induced protein, and is validated for WB, IHC-P, and ELISA. It cross-reacts with human, mouse, and rat, making it ideal for studying protein function in various signaling contexts.

[Learn More](#)

Parameter	Value
Product Name	FAM117A Rabbit pAb
Species	Human
Gene ID (human)	81558
Gene ID	81558
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 110-240 of human FAM117A (NP_110429.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9C073
Protein Weight	48kDa
Shipping	Ice bag