

Anti-Fam120A Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nzb2

Item Number: CM0021207

Introduction

Rabbit polyclonal antibody targeting human FAM120A (OSSA/C9orf10), suitable for WB and ELISA. Detects a ~122kDa membrane-associated scaffold protein involved in oxidative stress signaling, RNA metabolism, and lipogenic gene regulation.

[Learn More](#)

Parameter	Value
Product Name	FAM120A Rabbit pAb
Remarks/Alias	OSSA; C9orf10; HBVPTAP; FAM120A
Species	Human
GeneID	23196
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human FAM120A (NP_055427.2))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NZB2
Protein Weight	122kDa
Shipping	Ice bag