

Anti-Fam127A/Rtl8C Polyclonal Antibody For Wb And Elisa - A6Zki3

Item Number: CM0017862

Introduction

Rabbit polyclonal antibody against human FAM127A/RTL8C, validated for WB and ELISA, cross-reacts with mouse and rat. Ideal for retrotransposon Gag-like protein detection studies.

[Learn More](#)

Parameter	Value
Product Name	FAM127A Rabbit pAb
Remarks/Alias	CXX1; Mar8; MAR8C; Mart8; SIRH5; MART8C; FAM127A
Species	Human
Gene ID (Human)	8933
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-113 of human FAM127A (NP_001071639.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	A6ZKI3
Protein Weight	13kDa
Shipping	Ice bag