

Anti-Fam127B Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Bwd3

Item Number: CM0022414

Introduction

Rabbit polyclonal antibody against FAM127B (RTL8A/CXX1B/MAR8A), raised with recombinant full-length protein. Suitable for WB, IF/ICC, ELISA in human and mouse samples.

[Learn More](#)

Parameter	Value
Product Name	FAM127B Rabbit pAb
Remarks/Alias	CXX1b; MAR8A; SIRH6; FAM127B
Species	Human
GeneID (Human)	26071
GeneID	26071
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BWD3
Protein Weight	13kDa
Shipping	Ice bag