


Anti-Fam129B Rabbit Polyclonal Antibody For Wb, Elisa - Q96Ta1

Item Number: CM0021630



Introduction

FAM129B Rabbit Polyclonal Antibody validated for WB and ELISA, targeting human FAM129B (Niban 2). May play roles in apoptosis suppression and melanoma cell invasion. Suitable for cancer research.

[Learn More](#)

Parameter	Specification
Product Name	FAM129B Rabbit pAb
Remarks/Alias	OC58; MEG-3; C9orf88; FAM129B; MINERVA; bA356B19.6
Species	Human
GeneID (Human)	64855
GeneID	64855
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 447-746 of human FAM129B (NP_073744.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96TA1
Protein Weight	84kDa
Shipping	Ice bag