

Anti-Fam13A Polyclonal Antibody For Wb, Elisa - 094988

Item Number: CM0012948

Introduction

Rabbit polyclonal antibody against FAM13A, validated for WB and ELISA. Detects human FAM13A, a potential GTPase-activating protein with a predicted molecular weight of 117 kDa. Ideal for research on alternative splicing and phosphorylation.

[Learn More](#)

Parameter	Value
Product Name	FAM13A Rabbit pAb
Remarks/Alias	FAM13A1; ARHGAP48; FAM13A
Species	Human
GeneID (Human)	10144
GeneID	10144
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O94988
Protein Weight	117kDa
Shipping	Ice bag