

Anti-Fam155B Rabbit Polyclonal Antibody For Wb, Elisa - O75949

Item Number: CM0018740

Introduction

Rabbit polyclonal antibody targeting human FAM155B (NALF2/TED). Validated for Western blotting and ELISA; cross-reacts with mouse and rat. Supplied for ion channel and neuroscience research applications.

[Learn More](#)

Parameter	Specification
Product Name	FAM155B Rabbit pAb
Remarks/Alias	TED; TMEM28; CXorf63; FAM155B; bB57D9.1
Species	Human
GeneID (Human)	27112
GeneID	27112
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 400-472 of human FAM155B (NP_056501.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O75949
Protein Weight	52kDa
Shipping	Ice bag