

Anti-Fam160B2 Polyclonal Antibody For Wb, Elisa - Q86V87

Item Number: CM0026757

Introduction

Rabbit polyclonal antibody targeting human FAM160B2 (FHIP2B/RAI16), suitable for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Useful for investigating MAPK/ERK and TGFβ signaling pathways.

[Learn More](#)

Parameter	Description
Product Name	FAM160B2 Rabbit pAb
Remarks/Alias	FAI16; RAI16; FAM160B2; RAM160B2
Species	Human
GeneID	64760
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 420-520 of human FAM160B2 (NP_073586.5).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86V87
Protein Weight	82kDa
Shipping	Ice bag