

Anti-Fam173B Polyclonal Antibody For Wb And Elisa - Q6P4H8

Item Number: CM0027342

Introduction

Rabbit polyclonal antibody against human FAM173B (ATP synthase subunit C lysine N-methyltransferase). Validated for WB and ELISA, cross-reacts with human, mouse, and rat. Ideal for research on mitochondrial methylation and chronic pain.

[Learn More](#)

Parameter	Value
Product Name	FAM173B Rabbit pAb
Remarks/Alias	JS-2; FAM173B; hFAM173B
Species	Human
Gene ID (Human)	134145
Gene ID	134145
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-233 of human FAM173B (NP_954584.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6P4H8
Protein Weight	26kDa
Shipping	Ice bag