


Anti-Fam208A Rabbit Polyclonal Antibody For Wb, Elisa - Q9Uk61

Item Number: CM0012949



Introduction

Anti-FAM208A rabbit polyclonal antibody validated for WB and ELISA. Targets human TASOR, a key component of HUSH complex involved in epigenetic silencing of retrotransposons. Ideal for chromatin and gene regulation studies.

[Learn More](#)

Parameter	Value
Product Name	FAM208A Rabbit pAb
Remarks/Alias	RAP140; TASOR1; se89-1; C3orf63; FAM208A
Species	Human
GeneID (Human)	23272
GeneID	23272
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 700-1000 of human FAM208A (NP_001106207.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UK61
Protein Weight	189kDa
Shipping	Ice bag