


Anti-Fam20A Monoclonal Antibody For Wb - Q96Mk3

Item Number: CM0006726



Introduction

FAM20A Rabbit monoclonal antibody for Western blot (WB) and ELISA. Recognizes human pseudokinase FAM20A, a biomineralization regulator. Ideal for dental development and mineralization studies.

[Learn More](#)

Parameter	Value
Product Name	FAM20A Rabbit mAb
Remarks/Alias	AI1G; AIGFS; FP2747; FAM20A
Species	Human
GeneID (Human)	54757
GeneID	54757
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 442-541 of human FAM20A (NP_060035.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96MK3
Protein Weight	61 kDa
Shipping	Ice bag