


Anti-Fam20A Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q96Mk3

Item Number: CM0027641



Introduction

Rabbit polyclonal antibody targeting human FAM20A pseudokinase, validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Ideal for biomineralization and enamel research.

[Learn More](#)

Core Product Specifications	Parameter Details
Product Name	FAM20A Rabbit pAb
Remarks/Alias	AI1G; AIGFS; FP2747; FAM20A
Species	Human
GeneID (Human)	54757
GeneID	54757
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 140-240 of human FAM20A (NP_060035.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96MK3
Protein Weight	61kDa
Shipping	Ice bag