

Anti-Fam221A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8C790

Item Number: CM0023131

Introduction

Fam221a Rabbit pAb is a polyclonal antibody raised in rabbit against mouse FAM221A. Validated for WB, IF/ICC, ELISA, and cross-reactive with human and rat. Suitable for protein detection and localization studies.

[Learn More](#)

Parameter	Value
Product Name	Fam221a Rabbit pAb
Remarks/Alias	D330028D13Rik; Fam221a
Species	Mouse
Gene ID (Human)	340277
Gene ID	231946
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-292 of mouse Fam221a (NP_766315.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8C790
Protein Weight	33kDa
Shipping	Ice bag