

Anti-Fam3A Monoclonal Antibody For Wb, Ihc-P, Elisa - P98173

Item Number: CM0006475

Introduction

Recombinant rabbit monoclonal antibody targeting human FAM3A (P98173). Validated for WB, IHC-P, and ELISA; cross-reacts with mouse and rat. Ideal for defensin and fungal immunity studies.

[Learn More](#)

Parameter	Value
Product Name	FAM3A Rabbit mAb
Remarks/Alias	DLD; 2.19; XAP-7; DXS560S; FAM3A
Species	Human
Gene ID (Human)	60343
Gene ID	60343
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 34-230 of human FAM3A(NP_068578.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P98173
Protein Weight	25kDa
Shipping	Ice bag