

Anti-Fam3A Polyclonal Antibody For Wb, Ihc-P, Elisa - P98173

Item Number: CM0023700

Introduction

Rabbit polyclonal antibody targeting human FAM3A (aa 34-230), validated for WB, IHC-P, and ELISA. Recognizes a 25 kDa secreted cytokine-like protein with potential defensin activity.

[Learn More](#)

Parameter	Specification
Product Name	FAM3A Rabbit pAb
Remarks/Alias	DLD; 2.19; XAP-7; DXS560S; FAM3A
Species	Human
GeneID (Human)	60343
GeneID	60343
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 34-230 of human FAM3A (NP_068578.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P98173
Protein Weight	25kDa
Shipping	Ice bag