

Anti-Fam3B Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P58499

Item Number: CM0017871

Introduction

Rabbit polyclonal antibody against human FAM3B (PANDER), cross-reactive with mouse and rat. Validated for WB, IHC-P, IF/ICC, and ELISA. Used to study apoptosis in alpha and beta cells and pancreatic function.

[Learn More](#)

Parameter	Value
Product Name	FAM3B Rabbit pAb
Remarks/Alias	2-21; ORF9; PANDER; PRED44; C21orf11; C21orf76; FAM3B
Species	Human
Gene ID (Human)	54097
Gene ID	54097
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P58499
Protein Weight	26kDa
Shipping	Ice bag