

Anti-Fam3D Polyclonal Antibody For Wb, Elisa - Q96Bq1

Item Number: CM0020923

Introduction

CamelBio offers high-purity Anti-FAM3D Rabbit Polyclonal Antibody (Q96BQ1), validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for studying secreted glycoprotein FAM3D's role in placental biology and lectin-related pathways.

[Learn More](#)

Parameter	Value
Product Name	FAM3D Rabbit pAb
Remarks/Alias	EF7; OIT1; FAM3D
Species	Human
GeneID	131177
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-204 of human FAM3D (NP_620160.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q96BQ1
Protein Weight	25kDa
Shipping	Ice bag