

Anti-Fam3D Rabbit Polyclonal Antibody (Pab) For Wb, Ihc-P, If/Icc, Elisa - Q96Bq1

Item Number: CM0022955

Introduction

Rabbit polyclonal antibody against human FAM3D (EF7/OIT1) for Western blot, immunohistochemistry, immunofluorescence, and ELISA. Validated in human, mouse, and rat samples. Predicted molecular weight 25 kDa.

[Learn More](#)

Parameter	Value
Product Name	FAM3D Rabbit pAb
Remarks/Alias	EF7; OIT1; FAM3D
Species	Human
Gene ID (Human)	131177
Gene ID	131177
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96BQ1
Protein Weight	25kDa
Shipping	Ice bag