

Anti-Fam46A Polyclonal Antibody For Wb, Elisa - Q96Ip4

Item Number: CM0016513



Introduction

FAM46A rabbit polyclonal antibody validated for Western blot and ELISA. Detects human and mouse terminal nucleotidyltransferase 5A, a cytoplasmic poly(A) polymerase involved in bone mineralization and osteoblast function.

[Learn More](#)

Parameter	Value
Product Name	FAM46A Rabbit pAb
Remarks/Alias	OI18; XTP11; FAM46A; C6orf37
Species	Human
GeneID (human)	55603
GeneID	55603
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-400 of human FAM46A (NP_060103.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96IP4
Protein Weight	50kDa
Shipping	Ice bag