

Anti-Fam50B Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Wtj8

Item Number: CM0022146

Introduction

Rabbit polyclonal antibody against mouse FAM50B, suitable for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Useful for testis expression studies and acetylation research. Immunogen: recombinant fusion protein aa 70-150.

[Learn More](#)

Parameter	Value
Product Name	Fam50b Rabbit pAb
Remarks/Alias	X5L; D0H6S2654E; XAP-5-like; Fam50b
Species	Mouse
Gene ID (Human)	26240
Gene ID (Mouse)	108161
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 70-150 of mouse Fam50b (NP_620085.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9WTJ8
Protein Weight	40kDa
Shipping	Ice bag