

Anti-Fam65B Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q9Y4F9

Item Number: CM0026658

Introduction

Rabbit polyclonal antibody against human FAM65B (RIPOR2), validated for WB, IF-P, ELISA. Cross-reacts with mouse and rat. Ideal for studying cell polarization, myogenesis, and deafness-related pathways.

[Learn More](#)

Parameter	Value
Product Name	FAM65B Rabbit pAb
Remarks/Alias	PL48; DFNA21; DIFF40; DIFF48; FAM65B; MYONAP; C6orf32; DFNB104
Species	Human
GeneID (Human)	9750
GeneID	9750
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human FAM65B (NP_056948.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4F9
Protein Weight	119kDa
Shipping	Ice bag