

Anti-Fam96B Polyclonal Antibody For Wb, Elisa - Q9Y3D0

Item Number: CM0026458

Introduction

Rabbit polyclonal antibody against human FAM96B (CIA2B/MIP18) validated for Western blot and ELISA, also reactive with mouse. Suitable for iron-sulfur cluster assembly and chromosome segregation research.

[Learn More](#)

Parameter	Value
Product Name	FAM96B Rabbit pAb
Remarks/Alias	CIA2B; MIP18; FAM96B; CGI-128
Species	Human
GeneID (Human)	51647
GeneID	51647
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-163 of human FAM96B (NP_057146.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y3D0
Protein Weight	18kDa
Shipping	Ice bag