

## Anti-Fam98A Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Nca5

Item Number: CM0022204

### Introduction

Rabbit polyclonal antibody targeting human FAM98A, validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat samples. Ideal for studying PRMT1-induced arginine methylation and skeletal homeostasis.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | FAM98A Rabbit pAb                                                                                                 |
| Remarks/Alias    | FAM98A                                                                                                            |
| Species          | Human                                                                                                             |
| GeneID (Human)   | 25940                                                                                                             |
| GeneID           | 25940                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human FAM98A (NP_056290.3) |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, IHC-P, ELISA                                                                                                  |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | Q8NCA5                                                                                                            |
| Protein Weight   | 55kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |