

Anti-Fam98B Rabbit Polyclonal Antibody For Wb And Elisa - Q52Lj0

Item Number: CM0023789

Introduction

Rabbit polyclonal antibody against human FAM98B (tRNA-splicing ligase complex subunit) for WB and ELISA; cross-reacts with mouse and rat. Suitable for studying tRNA processing and colorectal cancer research.

[Learn More](#)

Parameter	Value
Product Name	FAM98B Rabbit pAb
Remarks/Alias	FAM98B
Species	Human
GeneID (Human)	283742
GeneID	283742
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-225 of human FAM98B (NP_775882.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q52Lj0
Protein Weight	46kDa
Shipping	Ice bag