


Anti-Fanca Monoclonal Antibody For Wb, Elisa - O15360

Item Number: CM0009221



Introduction

Rabbit monoclonal antibody targeting human FANCA protein (163 kDa), validated for Western blot and ELISA. Ideal for Fanconi anemia research, interstrand cross-link repair, and chromosome stability studies.

[Learn More](#)

Parameter	Value
Product Name	FANCA Rabbit mAb
Remarks/Alias	FA; FA1; FAA; FAH; FA-H; FACA; FANCH; FANCA
Species	Human
Gene ID (Human)	2175
Gene ID	2175
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2-104 of human FANCA (O15360).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O15360
Protein Weight	163kDa
Shipping	Ice bag