


Anti-Fanca Rabbit Polyclonal Antibody For Wb, Elisa - O15360

Item Number: CM0027425



Introduction

Rabbit polyclonal antibody targeting human FANCA. Validated for Western blot (WB) and ELISA. Recognizes the Fanconi anemia group A protein, crucial for DNA interstrand crosslink repair and chromosome stability. Suitable for Fanconi anemia and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	FANCA Rabbit pAb
Remarks/Alias	FA; FA1; FAA; FAH; FA-H; FACA; FANCH; FANCA
Species	Human
Gene ID (Human)	2175
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-275 of human FANCA (NP_000126.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O15360
Protein Weight	163kDa
Shipping	Ice bag