


Anti-Fancc Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q00597

Item Number: CM0021047



Introduction

High-quality rabbit polyclonal antibody against FANCC, validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat FANCC protein. Ideal for DNA damage repair and Fanconi anemia research.

[Learn More](#)

Parameter	Value
Product Name	FANCC Rabbit pAb
Remarks/Alias	FA3; FAC; FACC; FANCC
Species	Human
Gene ID (Human)	2176
Gene ID	2176
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 340-558 of human FANCC (NP_000127.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00597
Protein Weight	63kDa
Shipping	Ice bag