

Anti-Fancd2 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q9Bxw9

Item Number: CM0003860



Introduction

Rabbit monoclonal antibody against human FANCD2, suitable for WB, IP, and ELISA. Reacts with human, mouse, and rat. Ideal for studying DNA repair, Fanconi anemia, and chromosomal stability.

[Learn More](#)

Parameter	Specification
Product Name	FANCD2 Rabbit mAb
Remarks/Alias	FA4; FAD; FACD; FAD2; FA-D2; FANCD; FANCD2
Species	Human
GeneID (Human)	2177
GeneID	2177
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human FANCD2 (Q9BXW9).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BXW9
Protein Weight	164 kDa
Shipping	Ice bag