

Anti-Fancd2 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q9Bxw9

Item Number: CM0019028



Introduction

FANCD2 rabbit polyclonal antibody for WB, IF/ICC, IP, ELISA. Detects human and rat FANCD2, a key DNA repair protein. Suitable for Fanconi anemia, chromosome instability, and DNA damage response studies.

[Learn More](#)

Parameter	Value
Product Name	FANCD2 Rabbit pAb
Remarks/Alias	FA4; FAD; FACD; FAD2; FA-D2; FANCD; FANCD2
Species	Human
Gene ID	2177
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 540-720 of FANCD2 (NP_001018125.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9BXW9
Protein Weight	164kDa
Shipping	Ice bag