


Anti-Fancg Rabbit Polyclonal Antibody For Wb, Elisa - O15287

Item Number: CM0029164



Introduction

Rabbit polyclonal antibody against human FANCG (FAG/XRCC9). Validated for Western blot and ELISA, cross-reacts with human and mouse. Suitable for DNA repair and Fanconi anemia research.

[Learn More](#)

Parameter	Value
Product Name	FANCG Rabbit pAb
Remarks/Alias	FAG, XRCC9, FANCG
Species	Human
GeneID (Human)	2189
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O15287
Protein Weight	69kDa
Shipping	Ice bag