

Anti-Fanci Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nvi1

Item Number: CM0030172



Introduction

Rabbit polyclonal antibody targeting human FANCI, validated for Western blot (WB) and ELISA. Ideal for DNA damage repair, Fanconi anemia, and cell cycle checkpoint research.

[Learn More](#)

Parameter	Value
Product Name	FANCI Rabbit pAb
Remarks/Alias	KIAA1794; FANCI
Species	Human
Gene ID (Human)	55215
Gene ID	55215
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1061-1328 of human FANCI (NP_001106849.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NVI1
Protein Weight	149kDa
Shipping	Ice bag