


Anti-Fancl Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nw38

Item Number: CM0030124



Introduction

FANCL Rabbit Polyclonal Antibody validated for WB and ELISA, cross-reacts with mouse. Targets the E3 ubiquitin ligase involved in Fanconi anemia DNA repair. Ideal for studying DNA damage response and germ cell proliferation.

[Learn More](#)

Parameter	Value
Product Name	FANCL Rabbit pAb
Remarks/Alias	POG; PHF9; FAAP43; FANCL
Species	Human
GeneID (human)	55120
GeneID	55120
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9NW38
Protein Weight	43kDa
Shipping	Ice bag