


Anti-Fancm Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Iyd8

Item Number: CM0026906



Introduction

Rabbit polyclonal antibody targeting human FANCM, validated for Western blot, IF/ICC, and ELISA. Ideal for investigating Fanconi anemia pathway, DNA repair mechanisms, and germ cell biology.

[Learn More](#)

Parameter	Value
Product Name	FANCM Rabbit pAb
Remarks/Alias	POF15; SPGF28; FAAP250; KIAA1596; FANCM
Species	Human
Gene ID (Human)	57697
Gene ID	57697
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 390-660 of human FANCM (NP_065988.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q8IYD8
Protein Weight	232kDa
Shipping	Ice bag