

## Anti-Fap (Fibroblast Activation Protein-A) Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q12884

Item Number: CM0006339



### Introduction

Rabbit monoclonal antibody targeting human FAP (Q12884). Validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Essential tool for studying ECM remodeling, fibrosis, and tumor microenvironment.

[Learn More](#)

| Parameter        | Value                                                                                                     |
|------------------|-----------------------------------------------------------------------------------------------------------|
| Product Name     | Fibroblast activation protein- $\alpha$ (FAP) Rabbit mAb                                                  |
| Remarks/Alias    | FAP; DPPIV; FAPA; FAPalpha; SIMP; prolyl endopeptidase FAP; Fibroblast activation protein- $\alpha$ (FAP) |
| Species          | Human                                                                                                     |
| GeneID (Human)   | 2191                                                                                                      |
| GeneID           | 2191                                                                                                      |
| Immunogen        | Synthetic peptide corresponding to amino acids 411-510 of human FAP (NP_004451.2)                         |
| Source           | Rabbit                                                                                                    |
| Category         | Monoclonal Antibodies                                                                                     |
| Application      | WB, IF/ICC, ELISA                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                         |
| SWISS            | Q12884                                                                                                    |
| Protein Weight   | 26kDa/87kDa                                                                                               |
| Shipping         | Ice bag                                                                                                   |