

Anti-Fap Monoclonal Antibody For Flow Cytometry - Q12884

Item Number: CM0011707



Introduction

-conjugated rabbit monoclonal antibody targeting human Fibroblast activation protein alpha (FAP) for flow cytometry. Ideal for studying tissue remodeling, fibrosis, and tumor microenvironment. Validated for human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human FAP mAb
Remarks/Alias	FAP; DPPIV; FAPA; FAPalpha; SIMP; prolyl endopeptidase FAP
Species	Human
GeneID (Human)	2191
GeneID	2191
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q12884
Protein Weight	27kDa/87kDa
Shipping	Ice bag