

Anti-Fap Polyclonal Antibody For If/Icc, Elisa - Q12884

Item Number: CM0028573



Introduction

FAP Rabbit pAb is a human-specific polyclonal antibody for IF/ICC and ELISA, targeting fibroblast activation protein alpha (Q12884) involved in tumor growth and tissue remodeling. Ideal for cancer and fibrosis research.

[Learn More](#)

Parameter	Value
Product Name	FAP Rabbit pAb
Remarks/Alias	FAPA; SIMP; DPPIV; FAPalpha
Species	Human
GeneID (Human)	2191
GeneID	2191
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 481-580 of human FAP (NP_004451.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q12884
Protein Weight	88kDa
Shipping	Ice bag