

## Anti-Far1 Polyclonal Antibody For Wb, Elisa - Q8Wvx9

Item Number: CM0023639

### Introduction

Rabbit polyclonal antibody against human FAR1 (fatty acyl-CoA reductase 1), suitable for WB and ELISA. Targets a key enzyme in fatty alcohol, ether lipid, and wax monoester biosynthesis. Ideal for lipid metabolism and peroxisome research.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | FAR1 Rabbit pAb                                                                                                 |
| Remarks/Alias    | CSPSD; PFCRD; MLSTD2; SDR10E1; FAR1                                                                             |
| Species          | Human                                                                                                           |
| GeneID (human)   | 84188                                                                                                           |
| GeneID           | 84188                                                                                                           |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human FAR1 (NP_115604.1) |
| Source           | Rabbit                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                           |
| Application      | WB, ELISA                                                                                                       |
| Cross Reactivity | Human                                                                                                           |
| SWISS            | Q8WVX9                                                                                                          |
| Protein Weight   | 59kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                         |