

Anti-Farsb Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nsd9

Item Number: CM0017921



Introduction

Rabbit polyclonal antibody against human FARSF (phenylalanine-tRNA synthetase beta subunit), suitable for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studies on cytoplasmic aminoacylation and protein translation.

[Learn More](#)

Parameter	Value
Product Name	FARSF Rabbit pAb
Remarks/Alias	FRSB; RJSB; PheHB; PheRS; FARSLB; RILDBC; HSPC173; NEDBLLA; RILDBC1; FARSF
Species	Human
Gene ID (Human)	10056
Gene ID	10056
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 240-589 of human FARSF (NP_005678.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NSD9
Protein Weight	66kDa
Shipping	Ice bag