

Anti-Fas Monoclonal Antibody For Wb, Ihc-P, Elisa - P25445

Item Number: CM0003905



Introduction

Rabbit monoclonal antibody against human Fas (CD95, APO-1, TNFRSF6). Validated for WB, IHC-P, ELISA with cross-reactivity to rat. Ideal for apoptosis and immune tolerance studies.

[Learn More](#)

Parameter	Value
Product Name	Fas Rabbit mAb
Remarks/Alias	APT1; CD95; FAS1; APO-1; FASTM; ALPS1A; TNFRSF6; Fas
Species	Human
GeneID (Human)	355
GeneID	355
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human Fas (NP_000034.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P25445
Protein Weight	9-38 kDa
Shipping	Ice bag