

Anti-Fas Rabbit Polyclonal Antibody For Wb, Elisa - P25445

Item Number: CM0017590



Introduction

Rabbit polyclonal antibody targeting human Fas (CD95/APO-1), suitable for Western blot and ELISA. Cross-reacts with mouse and rat, making it ideal for apoptosis and immune regulation studies.

[Learn More](#)

Parameter	Value
Product Name	Fas Rabbit pAb
Remarks/Alias	APT1; CD95; FAS1; APO-1; FASTM; ALPS1A; TNFRSF6; Fas
Species	Human
GeneID (human)	355
GeneID	355
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 200-300 of human Fas (NP_000034.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P25445
Protein Weight	9-38 kDa
Shipping	Ice bag