


Anti-Fascin/Fscn1 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q16658

Item Number: CM0007787



Introduction

Rabbit monoclonal antibody against Fascin (FSCN1), validated for WB, IHC-P, ELISA in human, mouse, rat. Recognizes the 55 kDa actin-bundling protein involved in cell protrusions and migration. Ideal for cytoskeleton and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Fascin/FSCN1 Rabbit mAb
Remarks/Alias	HSN; SNL; p55; FAN1; Fascin/FSCN1
Species	Human
GeneID (human)	6624
GeneID	6624
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 249-362 of human Fascin/FSCN1 (Q16658).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16658
Protein Weight	55kDa
Shipping	Ice bag