

Anti-Fastkd1 Rabbit Polyclonal Antibody For Wb, Elisa - Q53R41

Item Number: CM0027762

Introduction

Rabbit polyclonal antibody against human FASTKD1, validated for Western blot and ELISA. Detects endogenous FASTKD1 in human, mouse, and rat samples. Suitable for mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	FASTKD1 Rabbit pAb
Remarks/Alias	FASTKD1
Species	Human
Gene ID (Human)	79675
Gene ID	79675
Immunogen	Recombinant fusion protein containing amino acids 638-847 of human FASTKD1 (NP_078898.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q53R41
Protein Weight	97kDa
Shipping	Ice bag