

Anti-Fastkd2 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Nyy8

Item Number: CM0030618

Introduction

High-quality rabbit polyclonal antibody targeting FASTKD2, validated for WB, IHC-P, and ELISA. Detects human and mouse FASTKD2 protein (81 kDa). Crucial for mitochondrial large ribosomal subunit assembly studies.

[Learn More](#)

Parameter	Value
Product Name	FASTKD2 Rabbit pAb
Remarks/Alias	COXPD44; KIAA0971
Species	Human
GeneID (Human)	22868
GeneID	22868
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 51-150 of human FASTKD2 (NP_055744.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NYY8
Protein Weight	81kDa
Shipping	Ice bag