


Anti-Fat10/Ubd Monoclonal Antibody For Wb And Elisa - O15205

Item Number: CM0007959



Introduction

Recombinant rabbit monoclonal antibody targeting human FAT10/UBD (Ubiquitin D). Validated for Western blot and ELISA, with cross-reactivity to mouse. Suitable for studies on ubiquitin-like protein degradation, NF-κB signaling, and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	FAT10/UBD Rabbit mAb
Remarks/Alias	FAT10; UBD-3; GABBR1; FAT10/UBD
Species	Human
GeneID	10537
Immunogen	Synthetic peptide corresponding to amino acids 70-150 of human FAT10/UBD (O15205)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O15205
Protein Weight	18 kDa
Shipping	Ice bag