


Anti-Fat10/Ubd Polyclonal Antibody For If/Icc, Elisa - O15205

Item Number: CM0017857



Introduction

Rabbit polyclonal antibody targeting human ubiquitin D (FAT10/UBD), validated for IF/ICC and ELISA with cross-reactivity to mouse and rat. Ideal for studying ubiquitin-like modification, proteasomal degradation, and innate immunity pathways.

[Learn More](#)

Parameter	Value
Product Name	FAT10/UBD Rabbit pAb
Remarks/Alias	FAT10; UBD-3; GABBR1; FAT10/UBD
Species	Human
Gene ID (Human)	10537
Gene ID	10537
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-165 of human FAT10/UBD (NP_006389.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15205
Protein Weight	18kDa
Shipping	Ice bag